



ICON INSTRUMENTS COMPANY

(Manufacturer of Laboratory Equipment)

IIC-118



MAGNETIC STIRRER WITH HOT PLATE

- Better low speed stirring with even small liquid volume.
- Better quality motor for higher Torque at low speeds.
- Stirring accomplished by small Teflon coated magnetic paddle.
- Paddle size as per Magnetic stirrer cap. Requirement.
- Accurate speed control & Temperature control.
- S. Steel 304 grade Top housing.
- Better quality Kanthal heaters all fire proof wiring connection better durability.
- Designed to reduce the fumes entrance in the unit.

Applications :

For Medical, Research Institute, Bio-Tech Pharmaceutical Agriculture, Food Industries

Model	IIC 118	IIC 118 A	IIC 118 B
Capacity (L)	1 Liter	2 Liters	5 Liters
Heater Wattage "W"	150W	300W	500W
Dimension (Apox.)	160 x 200 x 160	160 x 200 x 160	175 x 225 x 185

Technical Parameters :

Control Panel : Speed Regulator, Temp. Indicator & Mains
Material : CRC/MS sheet duly powder coated
Optional : All S. Steel body having S. Steel Top Magnetic Stirrer without hot plate having different dimension
Duty Cycle : Continuous
Power Consumption : 220/230V AC 50 HZ single Phase

**DIGITAL MODEL ALSO
AVAILABLE**

Works : Z-67, Okhla Industrial Area, Phase-II, New Delhi - 110020

Phone : +91-11-26383676 • E-mail : iconinstruments.india@gmail.com, sales@iconinstruments.co.in

Website: www.iconinstruments.co.in