



ICON INSTRUMENTS COMPANY

IIC-313

TABLET MAKING MACHINE

It is a small production machine, Single punch for making tablet upto 1000mg as per size of the die & punch 6,9 and 12mm.

Automatic compression & ejection device with each revolution of the wheel drive, automatic feed and removal of tablets. Built on a heavy cast body with eccentric pressure system, stainless steel hopper provided connected by a rubber tube, easy setting for thickness and hardness of tablets with nuts & locknuts and set of two spanners. Die and set of punches are of Chromium steel, fitted on Stand.



Features:

- Easily rotating wheel for making tablets.
- Compressing lever for adjustable hardness of tablets.
- Adjustable lower punch for adjusting the thickness.
- Stainless steel die and punch size 6, 9, 12mm.
- Electrical model can be provided in same with ½h.p motor.

Works : Z-67, Okhla Industrial Area, Phase-II, New Delhi - 110020

Phone : +91-11-26383676 • E-mail : iconinstruments.india@gmail.com, sales@iconinstruments.co.in

Website: www.iconinstruments.co.in